

Specific Test Parameters											
General		LIH	ISE	Range							
Test Name:	MALB		▼	◀	▶	Type:	Urine	▼	Operation	Yes	▼
Sample:	Volume	15	uL	Dilution	10	uL	Pre-Dilution Rate:	1			
Reagents:	R1 Volume	75	uL	Dilution	75	uL	Min OD	Max OD			
	R2 Volume	60	uL	Dilution	10	uL	L	H			
Wavelength:	Pri	340	▼	Sec.	800	▼	Reagent OD limit:				
Method:	END		▼				First L	-2.0	First H	2.5	
Reaction Slope:	+		▼				Last L	-2.0	Last H	2.5	
Measuring Point 1:	First	0		Last	27		Dynamic Range:				
Measuring Point 2:	First	0		Last	10		L	1.6	H	500	
Linearity:			%				Correlation Factor:				
No Lag Time:			▼				A	1.0000	B	0.0000	
							Onboard Stability Period:	61			

Specific Test Parameters										
General		LIH	ISE	Range						
Test Name:	MALB		▼	◀	▶	Type:	Urine	▼		
Value/flag	#		▼	Level L :	#		Level H :	#		
Normal Ranges:		Age L		Age H						
<input type="checkbox"/>	Sex	Year	Month	Year	Month	L	H			
<input type="checkbox"/>	1.	#	#	#	#	#	#	#	#	#
<input type="checkbox"/>	2.	#	#	#	#	#	#	#	#	#
<input type="checkbox"/>	3.	#	#	#	#	#	#	#	#	#
<input type="checkbox"/>	4.	#	#	#	#	#	#	#	#	#
<input type="checkbox"/>	5.	#	#	#	#	#	#	#	#	#
<input type="checkbox"/>	6.	#	#	#	#	#	#	#	#	#
<input type="checkbox"/>	7. None Selected									
<input type="checkbox"/>	8. Out of Range		L	H						
Panic Value:					Unit:	mg/L	Decimal places:	#		

Calibration Specific												
General		ISE										
Test Name:	MALB		▼	◀	▶	Type:	Urine	▼				
Calibration Type:	5AB		▼	FORMULA:	POLYGONAL		▼	Counts:	#	Process:	CONC	▼
Point 1:	Cal. No.	#	OD	CONC	+	FACTOR/OD-L	-2.0	FACTOR/OD-H	2.5			
Point 2:	#			+		-2.0		2.5				
Point 3:	#			+		-2.0		2.5				
Point 4:	#			+		-2.0		2.5				
Point 5:	#			+		-2.0		2.5				
Point 6:												
Point 7:												
1-Point Cal.Point				<input type="checkbox"/> With CONC-0								
MB Type Factor					Calibration Stability Period	14						

+ MEDICON Microalbumin Calibrator (pr. code. 1578-0852) or OLYMPUS ODR 3024 Microalbumin Calibrator (pr. code. 4478-0856)  
# User defined