## LYSOZYME CONTROL

For discrete analyzers

Product code: 1478-1925

Package: Level 1: 2 x 1ml  
Level 2: 2 x 1ml

Storage: 2-8°C

For *in vitro* use only

### INTENDED USE

Material for the quality control of the determination of Lysozyme with MEDICON reagents or BECKMAN COULTER OSR series reagents on MEDILYZER®, Diatron Pictus series, BECKMAN COULTER AU400/600/600-IVD/640/2700/5400 ADVIA 1200/1800/2400, DIMENSION® RxL /RxL HM /RxL Max/RxL Max HM/ Xpand/ Xpand Plus/Xpand Plus HM/ EXL/EXL LOCI Module, and every other automated biochemical analyzer that reference ranges have been determined and are mentioned in the appropriate leaflet. For *in vitro* use only.

### COMPOSITION

Lyophilized, recombinant human lysozyme in buffer.

### WARNINGS – PRECAUTIONS

- For *in vitro* use only.
- Avoid contact with skin, eyes and mucous membranes.

### PREPARATION

Remove the rubber stopper very carefully. As the material is packed under vacuum, sudden air influx may force material out of the vial. Reconstitute by adding exactly 1 ml of cold (2°-8°C) deionized or distilled water. Place the rubber stopper back on the vial, swirl gently 2-3 times and let it stand for 10 minutes. Mix gently to obtain a homogeneous mixture. Ideally use a hematology tubes mixer for 20 minutes. Avoid foaming. Use 30 minutes after water addition. Replace cap immediately after use and store at 2-8°C. Unsuitable storage, handling, or errors during the analytical procedure may give erroneous results.

### MATERIALS NEEDED BUT NOT PROVIDED WITH THE KIT

No reagents or materials needed, other than the automated biochemical analyzer and the Lysozyme reagents.

### STORAGE – STABILITY

Unopened, the control is stable up to the expiry date stated on the label at 2-8°C. After opening, the control is stable for 1 week, when stored tightly capped, at 2-8°C.

### TEST PROCESS

Refer to the user's manual of the analyzer for calibration and quality control process. Each laboratory should establish its own control frequency, however good laboratory practice suggests that controls be tested each date patient samples are tested and each time calibration is performed. The results obtained by any individual laboratory may vary from the given mean value but should fall within the corresponding acceptable ranges given in the enclosed table. If any trends or sudden shifts in values are detected, review all operating parameters. Each laboratory should establish guidelines for corrective actions to be taken if controls do not recover within the specified limits. Make sure that the LOT on the vial and on this insert is the same.

### WASTE DISPOSAL

Dispose of all material according to local guidelines for potentially infectious materials. Material safety data sheet is available by MEDICON on request.

### ASSIGNED VALUES – Lot specific

Refer to the table of assigned values. Acceptable range represents the maximum allowed deviation for one measurement only.

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