CLINICAL CHEMISTRY CONTROL LEV. 2

For discrete analyzers
 Product code: 1578-0902-12
 Packaging: 12 x 5 ml (lyoph.)
 Store at 2-8°C
 For *in vitro* use only

INTENDED USE

Material for the internal quality control of the results taken with MEDICON reagents on MEDILYZER®, Diatron Pictus series, BECKMAN COULTER AU400/600-IVD/640/2700/5400, ADVIA 1200/1800/2400, DIMENSION® Rxl/RxL HM/RxL Max/RxL Max HM/Xpand/Xpand Plus/Xpand Plus HM/EXL/EXL LOCI Module and all other analyzers the reference ranges for which have been determined and are recorded in the respective leaflet.

Contains the following analytes: ADA, ADA2, Alkaline Phosphatase, Alanine Aminotransferase, Aspartate Aminotransferase, Albumin, Amylase, Inorganic Phosphorus, Calcium, Creatinine, Lactate Dehydrogenase, Lactate, G-Glutamate Dehydrogenase, Creatine Kinase, Creatine Kinase-MB Activity (CKMB), Lithium, Magnesium, Total Protein, Acid and Prostatic Phosphatase, Urea, Uric Acid, Glucose, Iron, Triglycerides, UIBC/TIBC, Phospholipids, Direct and Total Bilirubin, Cholesterol, HDL, LDL, Cholinesterase, Potassium, Sodium and Chloride. For *in vitro* use only.

COMPOSITION

Lyophilized human serum with biological additives. The analyte concentrations / activities are LOT dependent. Contains preservative.

PREPARATION

Remove the rubber stopper very carefully. As the material is packed under vacuum, sudden air influx may force material out of the vial. Reconstitute by adding exactly 5ml of cold (2°-8°C) deionized or distilled water. Place the rubber stopper back on the vial, swirl gently 2-3 times and let it stand for 10 minutes. Mix gently to obtain a homogeneous mixture. Ideally use a hematology tubes mixer for 20 minutes. Avoid foaming. Use 30 minutes after water addition. Replace cap immediately after use and store at 2-8°C. Unsuitable storage, handling, or errors during the analytical procedure may give erroneous results.

For the quality control process refer to the user's manual of the analyzer.

Each laboratory should establish its own control frequency, however good laboratory practice suggests that controls should be tested each day patient samples are tested and each time calibration is performed. The results obtained by any individual laboratory may vary from the given mean value but should fall within the corresponding acceptable ranges given in the enclosed table. If any trends or sudden shifts in values are detected, review all operating parameters. Each laboratory should establish guidelines for corrective actions to be taken if controls do not recover within the specified limits. Make sure that the LOT on the vial and on this insert is the same.

STORAGE - STABILITY

Unopened, the content of vials is stable at 2-8°C up to the expiration date stated on the label. Reconstituted serum is stable for 8 hours at 25°C, 7 days at 4°C and 1 month when frozen once at -20°C.

Exceptions:

Bilirubin (if stored protected from light):

- 4 days at 4°C
- Do not store at 15-25°C
- Do not freeze

Acid Phosphatase: For Total and Prostatic - Acid Phosphatase, the material should be stabilized by adding 20 µL of Medicon's ACP stabilizer to 1 mL of the serum exactly 30 minutes after reconstitution. After stabilization Total and Prostatic Acid Phosphatase are stable for:

- 2 hours at 25°C
- 2 days at 4°C
- 1 month at -20°C when frozen once

Alkaline phosphatase:

Levels in the reconstituted serum will rise during storage. It is recommended that the reconstituted serum be allowed to stand for 1 hour at 25°C before measurement.

HDL- cholesterol / LDL- cholesterol:

Do not freeze

It is recommended to dispense the reconstituted control in aliquots and store at -20 °C for future use taking into consideration the above limitations.

MATERIALS NEEDED BUT NOT PROVIDED WITH THE KIT

- Automated biochemical analyzer
- Redistilled water
- Clinical chemistry reagents
- Common laboratory equipment

DETERIORATION

The control material should not be used:

- After the expiration date.
- When microbial growth is evident.
- After prolonged exposure to direct sunlight or high temperature.

WARNINGS - PRECAUTIONS

- For *in vitro* use only.
- Contains sodium azide (NaN₃) <0.1%. Avoid contact with eyes, skin and mucous membranes.
- Every donor used for the preparation of this control was tested and found negative for HbsAg, anti-HCV, anti-HIV 1 and 2 by FDA approved methods. However since no test method can offer complete assurance that infectious agents are absent, this product should be handled as any potentially infectious material.


ASSIGNED VALUES - Lot specific

Refer to the table of assigned values. Acceptable range represents the maximum allowed deviation for one measurement only.

WASTE DISPOSAL

This product contains sodium azide (NaN₃), which forms sensitive explosive compounds with copper or lead. Flush waste pipes with water after the disposal of undiluted reagent in order to avoid azide build up in the drain. Dispose of all waste material in accordance with local guidelines. Safety data sheet is available for professional use on request.

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