## PROTEIN STANDARD CALIBRATOR

For automated chemistry analyzers

Product code: 1578-1190

Package: 1x1 ml Level 1  
1x1 ml Level 2  
1x1 ml Level 3  
1x1 ml Level 4  
1x1 ml Level 5

Store at 2-8°C

For *in vitro* use

### INTENDED USE

MEDICON Protein Standard calibrator is a human plasma calibrator intended to be used for the quantitative determination of a1-antitrypsin, a2-macroglobulin, a1-acid glycoprotein, ceruloplasmin, haptoglobin, transferrin, IgA, IgG, IgM, C3, C4, prealbumin, kappa light chains, lambda light chains on biochemical analyzers with MEDICON reagents on MEDILYZER®, Diatron Pictus series, BECKMAN COULTER AU400/600/600-IVD/640/2700/5400, ADVIA 1200/1800/2400, DIMENSION® R<sub>x</sub>L /R<sub>x</sub>L HM /R<sub>x</sub>L Max/R<sub>x</sub>L Max HM/ Xpand/ Xpand HM/Xpand Plus/Xpand Plus HM/ EXL/EXL LOCI Module, MAXMAT and all other types of automated biochemistry analyzers. For *in vitro* diagnostic use only.

### COMPOSITION

Variable amounts of the following human proteins in a filtered, defibrinated, liquid human plasma base:

a1-antitrypsin, a2-macroglobulin, a1-acid glycoprotein, ceruloplasmin, haptoglobin, transferrin, IgA, IgG, IgM, C3, C4, kappa light chains, lambda light chains, Prealbumin stabilizers and preservatives. Actual concentrations are lot dependent. Please verify your calibrator lot before programming the analyzer.

### PREPARATION

Calibrators are ready to use. Swirl each vial gently prior use. Avoid foaming. Replace cap immediately after use and store at 2-8°C. Unsuitable storage, handling, or errors during the analytical procedure may give erroneous results.

### CALIBRATION STORAGE AND STABILITY

The calibrators are stable, unopened up to the stated expiry date when stored at 2-8°C. Once opened the calibrator is stable for 6 weeks stored at 2-8°C.

### MATERIALS REQUIRED AND NOT INCLUDED

No other materials are required than the chemistry analyzer, reagents and controls

### ASSIGNED VALUES – Lot specific

Immunturbidimetry (37 °C), traceable to IFCC CRM 470.

### DETERIORATION

The calibrator should not be used:

- After prolonged exposure to high temperatures.
- After deep freeze.

### PRECAUTIONS AND WARNINGS

- To avoid the possible built up of azide compounds, flush waste pipes with water after the disposal of calibrator. Dispose of all waste materials in accordance with local guidelines.
- Biological materials of human origin contained in this product were tested for Anti-HCV, HbsAg and Anti-HIV 1 and 2 using FDA approved methods and were found to be non-reactive. As there is no known test method that can offer complete assurance that products derived from human blood will not transmit infectious agents this product should be handled as a potentially infectious material.
- A Material Safety Data Sheet is available for this product by Medicon.

### WASTE DISPOSAL

This product contains sodium azide (NaN<sub>3</sub>), which forms sensitive explosive compounds with copper or lead. Flush waste pipes with water after the disposal of undiluted reagent in order to avoid azide build up in the drain pipes. Dispose of all waste material in accordance with local guidelines. Safety data sheet is available for professional use on request.

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