IMMUNOLOGY CONTROL 1

For discrete analyzers
 Product code: 1578-1195-04
 Package: 4 x 3 ml
 Store at 2-8°C
 For in vitro use

INTENDED USE

Material for the internal quality control of the performance of MEDICON special proteins reagents on MEDILYZER®, Diatron Pictus series, BECKMAN COULTER AU400/600/600-IVD/640/2700/5400, ADVIA 1200/1800/2400, DIMENSION® Rxl /Rxl HM /Rxl Max/Rxl Max HM/ Xpand/ Xpand HM/Xpand Plus/Xpand Plus HM/ EXL/EXL LOCI Module, and every other automated biochemical analyzer the reference ranges for which have been determined and are mentioned in the appropriate leaflet, or BECKMAN COULTER reagents on BECKMAN COULTER AU400/600/600-IVD/640/2700/5400 automated biochemistry analyzers. For in vitro use only.

COMPOSITION

Liquid, human serum based control material with additives and preservative. The concentrations/activities of the proteins are LOT depended.

WARNINGS-PRECAUTIONS

- For in vitro use only.
- Every donor used for the preparation of the control was tested and found negative for HbsAg, anti-HCV and anti-HIV 1 and 2 with FDA approved methods. As there is no known test method that can offer complete assurance that products derived from human blood will not transmit infectious agents this product should be handled as a potentially infectious material.
- Contains sodium azide (NaN₃) <0.1% as preservative. Avoid contact with skin, eyes and mucous membranes.

PREPARATION

Material ready to use. Swirl the vials gently to acquire homogenous mix. Avoid foaming. Replace cap immediately after use and store at 2-8°C. Unsuitable storage, handling, or errors during the analytical procedure may give erroneous results.

MATERIALS NEEDED BUT NOT PROVIDED

- Special proteins reagents
- Automated biochemical analyzer
- Common laboratory equipment.

STORAGE AND STABILITY

The control is stable, unopened up to the expiry date stated on the label when stored at 2-8°C. Once opened the control is stable for 1 month when stored tightly capped at 2-8°C.

DETERIORATION

The material should not be used

- After prolonged exposure to high temperature
- After the expiry date or when microbial growth is evident.

TEST PROCEDURE

Refer to the user's manual of the analyzer for calibration and quality control process. Each laboratory should establish its own control frequency, however good laboratory practice suggests that controls be tested each date patient samples are tested and each time calibration is performed. The results obtained by any individual laboratory may vary from the given mean value but should fall within the corresponding acceptable ranges given in the enclosed table. If any trends or sudden shifts in values are detected, review all operating parameters. Each laboratory should establish guidelines for corrective actions to be taken if controls do not recover within the specified limits. Make sure that the LOT on the vial and on this insert is the same.

WASTE DISPOSAL

Dispose of all material according to local guidelines for potentially infectious materials. Sodium azide forms explosive compounds with lead or copper. Flush waste pipes with plenty of water after disposal of undiluted control to avoid sodium azide build up in the drain. Material safety data sheet is available by MEDICON on request.

ASSIGNED VALUES

Refer to the table of assigned values. The range represents the maximum acceptable deviation for one measurement only

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