



# NEW IMPROVED PAD FOR DETERMINATION OF LEUCOCYTES

<b>BENEFIT:</b>	Significantly better colour difference between negative and 1st positive field
<b>ADVANTAGES:</b>	User friendly Improved accuracy
<b>PRINCIPLE</b>	Without changes

## NEW CONCENTRATION SCALE



## PREVIOUS CONCENTRATION SCALE



### Change: New sequence of zones

BLOOD ● BILI ● UBG ● KETONE ● GLUCOSE ● PROTEIN ● pH ● NITRITE ● **LEUCOCYTE** ● SG

### instead of

Blood ● **Leucocyte** ● Bili ● UBG ● Ketone ● Glucose ● Protein ● pH ● Nitrite ● SG



Note: This change refers only to PHAN<sup>®</sup> diagnostic strips with leucocytes! All other strips remained without any changes!



<b>NefroPHAN leuco</b>	10003352
<b>DekaPHAN leuco (pck. 50 strips)</b>	10003350
<b>DekaPHAN leuco (pck. 100 strips)</b>	10007386
<b>UndekaPHAN</b>	10003354

**PLIVA-Lachema Diagnostika s.r.o.**  
 Karásek 1, 621 33 Brno, Czech Republic  
 tel.: +420 541 127 244, fax: +420 541 127 627  
 e-mail: diagnostics@lachema.cz, www.lachema.com

zony 001 101207en

BLUE POINT OF QUALITY